

mis[®] | 4M[™] MATRIX[™]
BONE GRAFT CEMENT

4MATRIX


Characterization

An innovative synthetic bone grafting product, developed to simplify dental bone grafting procedures



4MATRIX

Product Composition



4MATRIX

Advantages



“All-in-one” syringe

containing the granulated powder and physiological saline for easy wetting



2-3 minutes

After 2-3 minutes setting period the cement is hardens to provide a 100% mechanically stable environment for facilitated membrane application (if required) and wound closure



Mixture

The fast and slow resorbing bone material mixture promotes optimal patient's bone regeneration. The BCS matrix completely resorb (within 4-10 weeks) and HA acts as a longer space maintainer (which resorb within months-years)

Main indications for using Biomaterials

Dehiscence




Extraction Sockets



Ridge Augmentation



Sinus Augmentation



4MATRIX

Packaging

One Package With
One Driver



"All in one" syringe
1cc Graft Cement

